## **Feet and Adaptors**












#### **Aluminium BK kit**







## SACH foot (pediatric)





Product Code	Side	Foot Length	Heel Height	Max. Load
0-01-FTT-L-12	Left	12cm	10mm	60kg
0-01-FTT-R-12	Right	12cm	10mm	60kg
0-01-FTT-L-13	Left	13cm	10mm	60kg
0-01-FTT-R-13	Right	13cm	10mm	60kg
0-01-FTT-L-14	Left	14cm	10mm	60kg
0-01-FTT-R-14	Right	14cm	10mm	60kg
0-01-FTT-L-15	Left	15cm	10mm	60kg
0-01-FTT-R-15	Right	15cm	10mm	60kg
0-01-FTT-L-16	Left	16cm	10mm	60kg
0-01-FTT-R-16	Right	16cm	10mm	60kg
0-01-FTT-L-17	Left	17cm	10mm	60kg
0-01-FTT-R-17	Right	17cm	10mm	60kg

#### **SACH foot (open toe)**

1-01-FTT







Product Code	Side	Foot Length	Heel Height	Max. Load
1-01-FTT-L-18	Left	18cm	10mm	60kg
1-01-FTT-R-18	Right	18cm	10mm	60kg
1-01-FTT-L-19	Left	19cm	10mm	60kg
1-01-FTT-R-19	Right	19cm	10mm	60kg
1-01-FTT-L-20	Left	20cm	10mm	90kg
1-01-FTT-R-20	Right	20cm	10mm	90kg
1-01-FTT-L-21	Left	21cm	10mm	90kg
1-01-FTT-R-21	Right	21cm	10mm	90kg
1-01-FTT-L-22	Left	22cm	10mm	95kg
1-01-FTT-R-22	Right	22cm	10mm	95kg
1-01-FTT-L-23	Left	23cm	10mm	95kg
1-01-FTT-R-23	Right	23cm	10mm	95kg
1-01-FTT-L-24	Left	24cm	10mm	95kg
1-01-FTT-R-24	Right	24cm	10mm	95kg
1-01-FTT-L-25	Left	25cm	10mm	95kg
1-01-FTT-R-25	Right	25cm	10mm	95kg
1-01-FTT-L-26	Left	26cm	10mm	100kg
1-01-FTT-R-26	Right	26cm	10mm	100kg
1-01-FTT-L-27	Left	27cm	10mm	100kg
1-01-FTT-R-27	Right	27cm	10mm	100kg
1-01-FTT-L-28	Left	28cm	10mm	100kg
1-01-FTT-R-28	Right	28cm	10mm	100kg
1-01-FTT-L-29	Left	29cm	10mm	100kg
1-01-FTT-R-29	Right	29cm	10mm	100kg
1-01-FTT-L-30	Left	30cm	10mm	100kg
1-01-FTT-R-30	Right	30cm	10mm	100kg



Product Code 1-01-FTT-01

1-01-FTT-01 connection plate for FTT (one size, suitable for 1-01-FTT, size 20cm-30cm)



Single axis foot 1-01-FT





Product Code	Side	Foot Length	Heel Height	Max. Load
1-01-FT-L-21	Left	21cm	10mm	80kg
1-01-FT-R-21	Right	21cm	10mm	80kg
1-01-FT-L-22	Left	22cm	10mm	80kg
1-01-FT-R-22	Right	22cm	10mm	80kg
1-01-FT-L-23	Left	23cm	10mm	80kg
1-01-FT-R-23	Right	23cm	10mm	80kg
1-01-FT-L-24	Left	24cm	10mm	80kg
1-01-FT-R-24	Right	24cm	10mm	80kg
1-01-FT-L-25	Left	25cm	10mm	80kg
1-01-FT-R-25	Right	25cm	10mm	80kg
1-01-FT-L-26	Left	26cm	10mm	90kg
1-01-FT-R-26	Right	26cm	10mm	90kg
1-01-FT-L-27	Left	27cm	10mm	90kg
1-01-FT-R-27	Right	27cm	10mm	90kg
1-01-FT-L-28	Left	28cm	10mm	90kg
1-01-FT-R-28	Right	28cm	10mm	90kg
1-01-FT-L-29	Left	29cm	10mm	90kg
1-01-FT-R-29	Right	29cm	10mm	90kg
1-01-FT-L-30	Left	30cm	10mm	90kg
1-01-FT-R-30	Right	30cm	10mm	90kg



1-01-FT-01 connection cap for Regal FT

Product Code	Size
1-01-FT-01-21	21cm
1-01-FT-01-22	22cm
1-01-FT-01-23	23cm
1-01-FT-01-24	24cm
1-01-FT-01-25	25cm
1-01-FT-01-26	26cm
1-01-FT-01-27	27cm
1-01-FT-01-28	28cm

DYNA foot 1-01-DYNA

New!



Product Code	Side	Foot length	Max. Load
1-01-DYNA-L-22	Left	22cm	75kg
1-01-DYNA-R-22	Right	22cm	75kg
1-01-DYNA-L-23	Left	23cm	75kg
1-01-DYNA-R-23	Right	23cm	75kg
1-01-DYNA-L-24	Left	24cm	75kg
1-01-DYNA-R-24	Right	24cm	75kg
1-01-DYNA-L-25	Left	25cm	80kg
1-01-DYNA-R-25	Right	25cm	80kg
1-01-DYNA-L-26	Left	26cm	80kg
1-01-DYNA-R-26	Right	26cm	80kg
1-01-DYNA-L-27	Left	27cm	80kg
1-01-DYNA-R-27	Right	27cm	80kg
1-01-DYNA-L-28	Left	28cm	80kg
1-01-DYNA-R-28	Right	28cm	80kg
1-01-DYNA-L-29	Left	29cm	80kg
1-01-DYNA-R-29	Right	29cm	80ka



Product Code	Description
1-01-DYNA-01	Fixed adaptor for DYNA foot

## Male socket adaptor (pediatric)

New!

0-02-A1SC 0-02-A1AC

Product Code	Weight	Material	Max. Load
0-02-A1SC	55g	SS	45kg
0-02-A1AC	33g	AL	45kg

#### Female socket adaptor (pediatric)

0-02-A5AC

0-02-A1SC/0-02-A1AC



New!

0-02-A5AC

Product Code	Weight	Material	Max. Load
0-02-A5AC	42g	AL	45kg

#### Female tube adaptor (pediatric)

0-02-A6AC





New!

0-02-A6AC

Product Code	Weight	Material	Max. Load
0-02-A6AC	36g	AL	45kg

#### **Angled tube adaptor (pediatric)**

0-02-A10AC



New!

Product Code	Weight	Material	Max. Load
0-02-A10AC	50g	AL	45kg



#### **Attachment block (pediatric)**

0-02-A22C



0-02-A22C

#### New!

Product Code	Weight	Material	Max. Load
0-02-A22C	98a	_	45ka

## **SACH** ankle adaptor

#### 1-02-2N1T/1-02-2N1S/1-02-2N1A



Product Code	Weight	Material	Max.Load
1-02-2N1T	92g	TI	125kg
1-02-2N1S	118g	SS	125kg
1-02-2N1A	98a	AL	100kg

## **SACH ankle adaptor (pediatric)**

#### 0-02-2N1SC/0-02-2N1TC



#### New!

Product Code	Weight	Material	Max.Load	
0-02-2N1TC	39g	TI	45kg	
0-02-2N1SC	62g	SS	45kg	



#### Single-axis ankle adaptor

#### 1-02-2ST/1-02-2SS/1-02-DB1







1-02-2SS



1-02-DB1 (bearing axis)

1-02-2SS-01G 1-02-2SS-01B

Product Code	Size	Weight	Material	Max.Load
1-02-2ST-S	22-25cm	196g	TI	125kg
1-02-2ST-L	26-30cm	208g	TI	125kg
1-02-2SS-S	22-25cm	273g	SS	125kg
1-02-2SS-L	26-30cm	296g	SS	125kg
1-02-DB1	-	173g	AL	100kg



1-02-2SS-01P

Heel cushion for Single-axis ankle adaptors

Product Code	Description	Hardness	Color
1-02-2SS-01G	Heel cushion	Hard	Green
1-02-2SS-01B	Heel cushion	Med	Black
1-02-2SS-01P	Heel cushion	Soft	Pink

#### **Tube adaptor**

#### 1-02-T3T/1-02-T3S



60mm

1-02-T3T/ 1-02-T3S



45mm

Product Code	Height	Weight	Material	Max. Load
1-02-T3T-45	45mm	48g	TI	100kg
1-02-T3T-60	60mm	56g	TI	100kg
1-02-T3T-75	75mm	67g	TI	100kg
1-02-T3S-45	45mm	84g	SS	100kg
1-02-T3S-60	60mm	99g	SS	100kg
1-02-T3S-75	75mm	120g	SS	100kg

## Tube with female pyramid adaptor

#### 1-02-T2A1T/1-02-T2A1S



1-02-T2A1S-210



1-02-T2A1T-410/ 1-02-T2A1S-410

Product Code	Length	Diameter	Weight	Tube material	Head material	Max. Load
1-02-T2A1T-210	210mm	Ø30mm	158g	AL	TI	100kg
1-02-T2A1S-210	210mm	Ø30mm	198g	AL	SS	100kg
1-02-T2A1T-410	410mm	Ø30mm	260g	AL	TI	100kg
1-02-T2A1S-410	410mm	Ø30mm	308g	AL	SS	100kg

# Tube with female pyramid adaptor 1-02-BO200/1-02-BO250/1-02-BO400 1-02-BO200 1-02-BO400

1-02-BO250

Product Code	Length	Diameter	Weight	Material	Max. Load
1-02-BO200	200mm	Ø30mm	150g	AL	100kg
1-02-BO250	250mm	Ø30mm	179g	AL	100kg
1-02-BO400	400mm	Ø30mm	336a	AL	100kg

#### **Tube with SACH ankle**

#### 1-02-BT200/1-02-BT400



Product Code	Length	Diameter	Weight	Material	Max. Load
1-02-BT200	200mm	Ø30mm	223g	AL	100kg
1-02-BT400	400mm	Ø30mm	336g	AL	100kg

## Tube (pediatric)

#### 0-02-BC150/0-02-BC300



Product Code	Length	Diameter	Weight	Material	Max. Load
0-02-BC150	150mm	Ø22mm	55g	AL	45kg
0-02-BC300	300mm	Ø22mm	125a	Al	45kg

#### Tube 1-02-B200/1-02-B400



Product Code	Length	Diameter	Weight	Material	Max. Load
1-02-B200	200mm	Ø30mm	118g	AL	100kg
1-02-B400	400mm	Ø30mm	298g	AL	100kg

Tube 1-02-C300/1-02-C355





1-02-C300

1-02-C355

Product Code	Length	Diameter	Weight	Material	Max. Load
1-02-C300	300mm	Ø30mm	115g	CF	100kg
1-02-C355	355mm	Ø30mm	133g	CF	100kg

#### Male socket adaptor

#### 1-02-A1T/1-02-A1S/1-02-A1A









1-02-A1T

1-02-A1S



Product Code	Weight	Material	Max. Load
1-02-A1T	53g	TI	125kg
1-02-A1S	120g	SS	125kg
1-02-A1A	59g	AL	100kg

1-02-A1T-01/	
1-02-A1S-01/ 1-02-A1A-01	
La	amination dummy

Product Code	Description
1-02-A1T-01	Lamination dummy
1-02-A1S-01	Lamination dummy
1-02-A1A-01	Lamination dummy

#### Socket adaptor

#### 1-02-TO1/1-02-TO2





1-02-TO1 1-02-T

Product Code	Weight	Material	Max. Load
1-02-TO1	160g	AL	100kg
1-02-TO2	170g	AL	100kg

#### Female socket adaptor

#### 1-02-A8T/1-02-A8S





1-02-A8T

1-02-A8S

 Product Code
 Weight
 Material
 Max. Load

 1-02-A8T
 57g
 TI
 100kg

 1-02-A8S
 96g
 SS
 100kg



## Female socket adaptor (Rotatable)

#### 1-02-A8ZT/1-02-A8ZS





1-02-A8ZS

1-02-A8ZT

Product Code	Weight	Material	Max. Load
1-02-A8ZT	65g	TI	100kg
1-02-A8ZS	117a	SS	100kg

#### **PU Foam block**

#### 1-02-A22



1-02-A22

Product Code	Weight	Material	Max. Load
1-02-A22	-	Foam	-

#### Wooden attachment block

#### 1-02-A22W



1-02-A22W

Product Code	Weight	Material	Max. Load
1-02-A22W	-	Wood	-

#### Male socket adaptor

#### 1-02-A2T/1-02-A2S/1-02-A2A





1-02-A2T

1-02-A2S



Product Code	Weight	Material	Max. Load
1-02-A2T	43g	TI	125kg
1-02-A2S	111g	SS	125kg
1-02-A2A	56g	AL	100kg



#### Male socket adaptor (Rotatable)

#### 1-02-A2ZT/1-02-A2ZS





1-02-A2ZT

1-02-A2ZS

Р	roduct Code	Weight	Material	Max. Load
	1-02-A2ZT	76g	TI	125kg
•	1-02-A2ZS	180g	SS	125kg

#### Female socket adaptor

#### 1-02-A5T/1-02-A5S/1-02-A5A





1-02-A5T

1-02-A5S



1-02-A5A

Product Code	Weight	Material	Max. Load
1-02-A5T	46g	TI	125kg
1-02-A5S	76g	SS	125kg
1-02-A5A	44g	AL	100kg

#### Female socket adaptor (Rotatable)

#### 1-02-A5ZT/1-02-A5ZS/1-02-TA1





1-02-A5ZT

1-02-A5ZS



Product Code	Weight	Material	Max. Load
1-02-A5ZT	90g	TI	125kg
1-02-A5ZS	142g	SS	125kg
1-02-TA1	102g	AL	100kg

#### Male socket adaptor

1-02-A12S



1-02-A12S

Product Code	Weight	Material	Max. Load
1-02-A12S	124g	SS	100kg



## Male socket adaptor (Rotatable)

1-02-A12ZT/1-02-A12ZS



Product Code	Weight	Material	Max. Load
1-02-A12ZT	145g	TI	125kg
1-02-A12ZS	168g	SS	125kg

## Female socket adaptor (Rotatable)

1-02-A11ZT/1-02-A11ZS



Product Code	Weight	Material	Max. Load
1-02-A11ZT	154g	TI	125kg
1-02-A11ZS	178g	SS	125kg

#### Anchor plate (KD)

0-00-4S41C=LA-01/0-00-4S40C=LA-01



0-00-4S41C=LA-01/ 0-00-4S40C=LA-01

#### New!

Product Code	Material	Description
0-00-4S41C=LA-01	SS	Anchor plate(KD)
0-00-4S40C=LA-01	SS	Anchor plate(KD)

#### Male tube adaptor

1-02-A7T/1-02-A7S







1-02-A7S

Product	Code	Weight	Material	Max. Load
1-02-	A7T	60g	TI	100kg
1-02-	47S	96g	SS	100kg



#### Female tube adaptor

#### 1-02-A6T/1-02-A6S/1-02-A6A





1-02-A6T

1-02-A6S



1-02-A6A

Product Code	Weight	Material	Max. Load
1-02-A6T	63g	TI	125kg
1-02-A6S	113g	SS	125kg
1-02-A6A	78a	AL	100kg

#### **Offset Tube Clamp**

#### 1-02-A6TOFF



1-02-A6TOFF

#### 1-02-A6TOFF **Offset Tube Clamp**

- The titanium 30mm offset tube clamp is used when extra alignment is needed
- It is also ideal when using a pin suspension lock, where prosthetic and pin alignment can be easily adjusted using the offset clamp
- This versatile clamp design also provides an offset position that can be rotated 360 degrees

Product Code	Weight	Material	Max. Load
1-02-A6TOFF	109g	TI	125kg

#### Angled tube adaptor

#### 1-02-A10T/1-02-A10S



1-02-A10T



1-02-A10S

Product Code	Weight	Material	Max. Load
1-02-A10T	73g	TI	125kg
1-02-A10S	125a	SS	125kg

#### **Rotation adaptor**

#### 1-02-A13T



1-02-A13T

#### 1-02-A13T Rotation adaptor

- Set Screws providing the lower connection
- With the incorporation of the rotation adapter above the knee
- The shin may be rotated relative to the socket with the knee flexed
- For rotation, press the black button then turn the knee

Product Code	Weight	Material	Max. Load
1-02-A13T	198g	TI	100kg



#### Sliding tube clamp

1-02-A14T/1-02-A14S



1-02-A14T/1-02-A14S Sliding tube clamp

- 30mm sliding tube clamp adaptor
- Permits additional sliding adjustment between the prosthetic socket and tube adaptor, even as the prosthetic is being worn
- The sliding adjustment allows for fine turning in any desired parallel position of the distal portion, independent og the angle of adjustment through the pyramid

Product Code	Weight	Material	Max. Load
1-02-A14T	180g	TI	100kg
1-02-A14S	313a	SS	100kg

## Double adaptor with adjusting of length

1-02-A15T/1-02-A15S



1-02-A15T/1-02-A15S

## 1-02-A15T/1-02-A15S Double adaptor with adjusting of length

- · Suitable for TT and TF amputee with long residual limb
- The height could be adjusted from 60mm to 75mm

Product Code	Weight	Height	Material	Max. Load
1-02-A15T	133g	60-75mm	TI	100kg
1-02-A15S	205g	60-75mm	SS	100kg

#### **Double adaptor**

1-02-A3T/1-02-A3S



\*\*





60mm

1-02-A3T/ 1-02-A3S

Product Code	Weight	Material	Max. Load
1-02-A3T-32	63g	TI	125kg
1-02-A3T-45	72g	TI	125kg
1-02-A3T-60	84g	TI	125kg
1-02-A3T-75	90g	TI	125kg
1-02-A3S-32	115g	SS	125kg
1-02-A3S-45	144g	SS	125kg
1-02-A3S-60	174g	SS	125kg
1-02-A3S-75	202a	SS	125ka



**Double adaptor** 

1-02-A4T/1-02-A4S







60mm

1-02-A4T/ 1-02-A4S

Product Code	Weight	Material	Max. Load
1-02-A4T-32	70g	TI	125kg
1-02-A4T-45	82g	TI	125kg
1-02-A4T-60	93g	TI	125kg
1-02-A4T-75	103g	TI	125kg
1-02-A4S-32	121g	SS	125kg
1-02-A4S-45	163g	SS	125kg
1-02-A4S-60	192g	SS	125kg
1-02-A4S-75	222g	SS	125kg

#### Single axis ankle joint (for wooden single axis ankle)

1-02-2SSW



New!

**Product Code** 

1-02-2SSW